

**WORK INSTRUCTION FOR KODAMA CATHETER's SPC ANALYSIS APPLICATION**

WI: 000-5-WI-134

Version: 1

Page: 1/3

Effective date: Jun 11th, 2012

**I. PURPOSE:**

This Work instruction provide the parameter which are must applied SPC control for Kodama catheter product.

**II. APPLICATION:**

This Work instruction concerns PRD, QAS Department

**III. REFERENCE DOCUMENTS:**

4-Pr-008:	SPC procedure
4-OP-345:	Kodama catheter
4-OP-614:	Hydrophilic coating

**IV. TERM DEFINITION****V. FLOWCHART**

UNCONTROL COPY IF PRINTOUT

Checked by : Dao Ngoc Trung

Date: 11 Jun 2012



Approved by: Nguyễn Hữu Phụng Vân

Date: 11 Jun 2012



Prepared by : Nguyen Thanh Trung

Date: 11 Jun 2012

# WORK INSTRUCTION FOR KODAMA CATHETER's SPC ANALYSIS APPLICATION

WI: 000-5-WI-134

Version: 1

Page: 2/3

## I. Scope:

This Working instruction apply for Kodama catheter product which is made in FOV.

## II. SPC control parameter:

NO	Location	Caterogy/Process	Parameter	Specification	Reference document	Remark
1	Medical cleaning room	Main cleaning room enviroment control	Particle	Size's 0.5 um .  Quantity < 5000 pcs/cfu.		
2			Temperature	20oC ~ 28oC.		
3			Humidity	10% ~ 90%.		
4		All heating processes	Temperature	Refer to process specifiati  in Operation procedure.	4-OP-345	
5		UV curing	UV power	Power ≥ 4 W/cm2		EFOS 100SS
6		Drive cable assembly	Laser power	0.0855 J ~ 0.1088J		
7	Hydrophilic coating room	Cleaning room enviroment control	Particle	Size's 0.5 um .  Quantity < 5000 pcs/cfu.		
8			Temperature	18oC ~ 25oC		
9			Humidity	Below 60%		
10		Hydrophilic coating	Viscosity	20% ~ 80%	4-OP-614	
11			Moisture	-		Waiting SVMI

**WORK INSTRUCTION FOR REWORK PRODUCT SVMI****WI: 000-5-WI-134****Version: 1****Page: 3/3****VI. REVISION HISTORY****REVISION HISTORY**

Date	Person	Version	Old content of description	New content		Reason of change
				New description	Change requester	
21-May-12	T.Trung	1		* Establish		

**UNCONTROL COPY IF PRINTOUT**